


Patient Care Guide

Get to know HelixBand


Wearing HelixBand



1. Secure your AFO

Insert your foot into your shoe with the HelixBand straps under your foot. The straps will be between your foot and the foot-plate of the AFO.

Secure and tighten the calf strap.



2. Tightening the strap

While seated, wrap the strap around your leg. Bend forward and pull tight. You will feel your arch being lifted.

Secure the Velcro hook to the loop area on the calf cuff.



3. Storage

When not using your AFO, it's helpful to keep your straps attached to the loop area and off the ground.

Other information

Cleaning information

Straps:

- Hand wash with mild detergent.
- Pat surface dry with towel and leave to air dry completely.
- Do not use bleach, fabric softener, chlorinated water, or tumble dry in a machine dryer.

Videos and accessories

Our website is a great resource to learn more about your Elevate product.

Find videos on how to properly wear your AFO and more. You can also buy extra liners, straps, and Elevate gear!

Learn more:
elevatemovement.com

Product issues?

If you have any issues with your product, please contact the clinician or clinic who performed your initial fitting to start a warranty claim with the Elevate team.

077-1031-A • Rev 01/23

More questions? We're here to help.
Reach out to us today or find more info online.

(888) 504-0230

Monday - Friday | 9am-8pm EST
www.elevatemovement.com

