

# Turtlebrace

## Fact sheet

### Turtlebrace for the wrist

Direct molding immobilization brace, for the wrist, made with low temp thermoplastic core.

#### Product

---

An “apply and forget” type brace made of stratification of multi-perforated low-temperature thermoplastic and a microfiber fabric.

The great flexibility of the components stabilizes the joint more easily, and with more precision, than a traditional brace or cast.

It can be reshaped many times, without losing its qualities, allowing adjustment according to the loss of swelling or morphological changes.



#### Characteristics

---

- Direct molding on the limb
- Antibacterial fabric
- Ambidextrous, it comes with Velcro or zipper closure
- Radiotransparent
- Submersible and auto-draining
- Remoldable without limits nor degradation
- 3.2 mm thickness that permits normal dressing
- Compostable eco-friendly plastic

#### Classification


---

Single patient class 1 medical device (FDA approved and CE certified)

#### Indications

---

Stable fractures, sprains, post-op immobilization, carpal tunnel, tendinitis, serial casting, spasticity control.




# Turtlebrace

## Models

---

	Pediatric	
	Velcro	Zipper
Medium (10.4 cm – 16.6 cm)	TBPPV-02	TBPPZ-02
Large (12.8 cm – 20.5 cm)	TBPPV-03	TBPPZ-03
	Adult	
	Velcro	Zipper
Small (15.2 cm – 24.3cm)	TBPAV-01	TBPAZ-01
Medium (17.6 cm – 28.2 cm)	TBPAV-02	TBPAZ -02
Large (20 cm – 32 cm)	TBPAV-03	TBPAZ -03



## Molding instruction

---

- 1- Heat the brace (dry heat or hot water), between 67 °C and 108 °C (152 °F and 225 °F) until it becomes soft and elastic. The plastic must fill doughy when pinch between two fingers.
  - Dry heat method  
Put the brace in either the Turtlebrace heating bag, a regular or convection oven. If you use a regular or convection oven, pre-heat them to 102 °C (215 °F) before heating the brace.
  - Hot water method  
Place the brace in an hydrocollator or a hot-water heating pan, between 67°C and 100°C (152°F and 212°F). If you use a hot water pan, make sure that the brace doesn't touch the bottom because the bottom temperature can exceed 108 °C (225 °F).
- 2- Once the brace has become soft and elastic, you can drape the brace on the body. Make sure that the temperature of the brace is not too hot for comfort or at risk of burning your patient. Small padding, about 3mm (1/8"), can be place at the bony apex and removed after the molding.
- 3- Place your patient in the desired position and wait for the brace to harden. To avoid rippling of the zipper, keep a tension on the top end of the zipper, or on both ends.

Once hardened, the brace will keep its shape until you heated back again.


## Recommendations

---

Molding should be done only by a health professional, or somebody trained in bracing, casting, or similar medical devices.

This is a single patient use, it cannot be transferred, even if it had been washed thoroughly.

Do not use a heat gun as it may burn the brace.



# Turtlebrace

## **Precautions**

---

Do not drape the brace if it is too hot to avoid skin burns or discomforts.

It is recommended to check the blood circulation often. If the brace becomes too tight, advise the client to unzip, loosen, or remove if possible, the brace and call his health professional.

It is recommended to check the skin often. If the show signs of maceration, irritation (redness), rashes, or other skin problem, advise the client to remove the brace (if possible) and immediately call his health professional.

Do not heat the brace over 108 °C (225 °F), because the fabric or/and the zipper could burn or melt.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **Information**

---

If you have any questions or in need of more information.

- Consult our website at [www.turtlebrace.com](http://www.turtlebrace.com) or

Write to: [info@turtlebrace.com](mailto:info@turtlebrace.com)