

A/K SHRINKERS

2



Patent pending.

COMPRESSOGRIP® A/K SHRINKERS


The new Compressogrip A/K for transfemoral applications provides the same advantages you have always enjoyed with the original Compressogrip Shrinker: clinically effective, easy to use and apply, lower cost than traditional shrinkers, fewer sizes fit more patients – easy and affordable to keep on hand for same day delivery.

The fabric used for these shrinkers is circular knit with increased power in the elastic, producing a heavier fabric than that used for B/K Compressogrip shrinkers. Multi-directional stretch for improved fitting qualities and a universal waist belt. Available in standard white Compressogrip and with X-Static® – The Silver Fiber®.

This fabric gives 25-30 mmHg compression at 50% stretch.

WHITE	X-STATIC	SIZE	FITS CIRCUMFERENCE		FLAT WIDTH
2H-OS05-XS	2H-0X05-XS	X-Small	12 – 17"	30 – 43 cm	5"
2H-OS06-SM	2H-0X06-SM	Small	15 – 20"	38 – 50 cm	6"
2H-OS07-MD	2H-0X07-MD	Medium	18 – 24"	45 – 60 cm	7.5"
2H-OS09-LG	2H-0X09-LG	Large	24 – 32"	60 – 80 cm	9"
6W-7322-30 — Auxiliary Suspension Laces, Sold 12 per package					

3-DIMENSIONAL A/K SHRINKER — Standard Sizes



Easy to adjust, hook and pile, elongate side opening waist band, adaptable for left or right amputation.

Optional auxiliary support for proximal edge at groin can be attached at any point on waist band or easily detached and removed.

Tapered design for gradient pressure; heavy or medium pressure available.

Superior 3-dimensional support and shaping of distal end; garment length from groin to distal end should be one to two inches shorter than actual measurement so that proper snugging of distal end will not result in roll at groin.

FLAT MEASURE		FITS CIRC.		GROIN to DISTAL LENGTH	60-65 mmHg	30-35 mmHg
TOE	TOP	TOE	TOP		HEAVY	MEDIUM
4"	6"	12"-14.5"	15"-17"	7"	2HPH6407	2HPM6407
				10"	2HPH6410	2HPM6410
				13"	2HPH6413	2HPM6413
4"	7"	17.5"-20.5"		7"	2HPH7407	2HPM7407
				10"	2HPH7410	2HPM7410
				13"	2HPH7413	2HPM7413
4"	8"	21"-23"		7"	2HPH8407	2HPM8407
				10"	2HPH8410	2HPM8410
				13"	2HPH8413	2HPM8413
5"	7"	15"-17.5"	17.5"-20.5"	7"	2HPH7507	2HPM7507
				10"	2HPH7510	2HPM7510
				13"	2HPH7513	2HPM7513
5"	8"	21"-23"		7"	2HPH8507	2HPM8507
				10"	2HPH8510	2HPM8510
				13"	2HPH8513	2HPM8513
5"	9"	23.5"-25.5"		7"	2HPH9507	2HPM9507
				10"	2HPH9510	2HPM9510
				13"	2HPH9513	2HPM9513
6"	8"	18"-22"	21"-23"	7"	2HPH8607	2HPM8607
				10"	2HPH8610	2HPM8610
				13"	2HPH8613	2HPM8613
6"	9"	23.5"-25.5"		7"	2HPH9607	2HPM9607
				10"	2HPH9610	2HPM9610
				13"	2HPH9613	2HPM9613
6"	10"	26"-29"		7"	2HPH0607	2HPM0607
				10"	2HPH0610	2HPM0610
				13"	2HPH0613	2HPM0613

CARE INSTRUCTIONS: Attach hook to pile before washing. Machine wash warm, no bleach; dry low.