

Modular Dynamic AFO

Care & Use

To reshape, spot heat the carbon fiber with a heat gun up to 250° F (121° C) (higher temperatures will cause delamination). **Do not heat or reshape the PEEK rods or the carbon fiber PEEK attachments.** Attach padding kit. Select a laced shoe with a 5/8 in. (16 mm) heel and removable insole. Make sure to place the foot plate under the insole or foot orthotic in the shoe. For best results, use a custom foot orthotic.

Warranty

1 year warranty*

Swedish Patent No. 531 069

U.S. Patent No. 7,766,851

Fillauer®

www.fillauer.com

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946
Fax 423.629.7936



© 2017 Fillauer LLC
PM295 07-11-12/01-31-17