



## F3 HERO PEDIATRIC AFO

Dynamic Feedback Anterior  
Ankle Foot Orthosis

OPTIMALLY REDESIGNED FOR  
PEDIATRIC USE



# F3 HERO AFO

By Thrive Orthopedics

484-442-0494 | [THRIVEORTHOPEDICS.COM](http://THRIVEORTHOPEDICS.COM)

## PEDIATRIC INDICATIONS

### AMBULATORY WEAKNESS ASSOCIATED WITH

Footdrop, idiopathic toe walker with no mid-foot collapse, gait deviations secondary to proprioceptive deficit, muscular dystrophy, spastic diplegia

## BENEFITS OF THE F3 HERO AFO

### FORM

Handcrafted and elegantly finished, striking the perfect balance between lightweight construction and durability.

### FIT

Made with the end-user in mind, Thrive's F3 Hero AFO comes in four convenient sizes with trimmable footplates for an ergonomic and intimate fit.

### FUNCTION

Strategically designed to enhance dynamic feedback throughout the gait cycle.

## F3 HERO AFO FITTING GUIDE

			
SKU	SIZE	HEIGHT	FOOTPLATE LENGTH
PASLT OR PASRT	SMALL	8 5/8"	6 1/4"
PAMLT OR PAMRT	MEDIUM	10 1/8"	7"
PALLT OR PALRT	LARGE	11 5/8"	7 7/8"
PAXLLT OR PAXLRT	EXTRA LARGE	12 3/8"	8 1/4"

FITTING RECOMMENDATIONS ARE BASED ON PREFABRICATED SIZES WITHOUT MODIFICATION

484-442-0494 | [THRIVEORTHOPEDICS.COM](https://www.thriveorthopedics.com)