



## **F3 AFO**

Dynamic Feedback Anterior  
Ankle Foot Orthosis

**OPTIMALLY DESIGNED FOR  
PATIENTS WITH DEXTERITY ISSUES**



484-442-0494 | [THRIVEORTHOPEDICS.COM](https://THRIVEORTHOPEDICS.COM)

## INDICATIONS

### AMBULATORY WEAKNESS ASSOCIATED WITH

- ✓ Foot Drop
- ✓ Foot Slap
- ✓ Multiple sclerosis
- ✓ Muscular dystrophy
- ✓ Traumatic brain injuries
- ✓ Charcotmarie-tooth
- ✓ Achilles Tendinosis
- ✓ Guillian-barre syndrome
- ✓ Post polio syndrome
- ✓ Stroke
- ✓ Cerebral palsy
- ✓ Myelomeningocele
- ✓ Mid-foot and ankle arthrosis

## BENEFITS OF THE F3 AFO

### FORM

Handcrafted and elegantly finished, striking the perfect balance between lightweight construction and durability.






### FIT

Made with the end-user in mind, Thrive's F3 AFO comes in five convenient sizes with trimmable footplates for an ergonomic and intimate fit.

### FUNCTION

Strategically designed to enhance dynamic feedback throughout the gait cycle.

## F3 AFO FITTING GUIDE

					
PART NUMBER	SIZE	US WOMEN'S SHOE SIZE	US MEN'S SHOE SIZE	FOOTPLATE LENGTH (inches)	BRACE HEIGHT (inches)
AXSLT OR AXSRT	XS	6-8.5	5-7	9"	13"
ASLT OR ASRT	S	9-10	7.5-9.5	9 <sup>3/4</sup> "	14"
AMLT OR AMRT	M	10.5-12	10-11	10 <sup>3/8</sup> "	15"
ALLT OR ALRT	L	12.5-14	11.5-13	11"	16"
AXLLT OR AXLRT	XL	14.5-16	13.5-15	11 <sup>3/4</sup> "	17"

FITTING RECOMMENDATIONS ARE BASED ON PREFABRICATED SIZES WITHOUT MODIFICATION

484-442-0494 | THRIVEORTHOPEDICS.COM