



INDICATIONS

AMBULATORY WEAKNESS ASSOCIATED WITH

- **∀** Foot Drop
- **ਓ** Foot Slap
- **⊗** Multiple sclerosis
- **⊘** Muscular dystrophy
- **⊗** Traumatic brain injuries
- *⊙* Charcotmarie-tooth
- **⊗** Achilles Tendinosis

- **⊘** Guillian-barre syndrome
- **⊘** Post polio syndrome
- **⊗** Stroke
- *⊙* Cerebral palsy
- ✓ Myelomeningocele
- ✓ Mid-foot and ankle arthrosis

BENEFITS OF THE F5 AFO

FORM

Handcrafted and elegantly finished, striking the perfect balance between lightweight construction and durability.

FIT

Made with the end-user in mind, Thrive's F5 AFO comes in five convenient sizes with trimmable footplates for an ergonomic and intimate fit.

FUNCTION

Strategically designed to enhance dynamic feedback throughout the gait cycle.

F5 AFO FITTING GUIDE

	. 1- 36		- 1.s. [3]	
	US WOMEN'S SHOE SIZE	US MEN'S SHOE SIZE	FOOTPLATE LENGTH (inches)	BRACE HEIGHT (inches)
xs	5.5 - 7.5	4.5 - 6.5	8 ^{7/8}	12
S	8 - 9.5	7 - 8.5	9 ^{5/8}	13
М	10 - 12	9 - 11	10 1/2	14
L	12.5 - 14.5	11.5 - 13.5	11 1/4	15
XL	15 +	14 - 16	12	16

FITTING RECOMMENDATIONS ARE BASED ON PREFABRICATED SIZES WITHOUT MODIFICATION