

## What providers and patients are saying about GlideWear®

"I recommend GlideWear® socks to at least five patients every day. Strategic friction reduction is a game changer for many patients with podiatric conditions that are not entirely alleviated by traditional treatments. Compliance is high since patients can tell right away that they are working."

*Troy J. Boffeli, DPM, FACFAS*

Department Chair and Residency Program Director  
HealthPartners Medical Group and Regions Hospital  
St. Paul, MN

"I have found the GlideWear® socks to be of excellent quality. I have a particular patient with bilateral foot ulcers, and am having him wear just the GlideWear socks in his regular shoes. He returns today saying he loved them, found them very comfortable, and is walking more than he would normally. Despite walking more, his wounds are much smaller as well!"

*Podiatrist from New Mexico*

"I have developed, from time to time, a dermatitis on top of my foot and around the ankle. The skin was sensitive not only during the day from my sock but also from the bed sheet. The heel and ankle sock relieved the symptoms and aided in a quick recovery of the area. If and when the dermatitis returns I'll grab my GlideWear® socks quickly!"

*Don H.*

"I have problems with calluses and I am diabetic with neuropathy. I would try anything to help my feet and what a difference these socks make. My feet feel so good, they are so comfortable!!! Love them!!"

*Sharon S.*



For more information about our products, visit

**[www.tamarackhti.com](http://www.tamarackhti.com)**

Email: [info@tamarackhti.com](mailto:info@tamarackhti.com)

Call: (763) 795-0057

Toll-Free: (866) 795-0057

Consumers may purchase GW-SX Socks via  
**Orthopedic Supplies of New Mexico.**

Practitioners or facilities may purchase through  
**Becker Orthopedic.**

**NOTICE:** GlideWear® is not intended for application directly against open wounds. Consult your doctor if you have any questions about the use of GlideWear® products with other prescribed therapies.



© 2023 Tamarack Habilitation Technologies, Inc.  
GlideWear is a registered trademark of MIPS, manufactured in the USA by Tamarack HTI.  
Patented US8646459 - US9241853 - GB2478206 - CA2738162 - Other Patents Pending. All Rights Reserved.

Tamarack Habilitation Technologies, Inc.  
1670-94th Lane NE, Blaine, MN 55449  
Tel: 866.795.0057  
E-mail: [info@tamarackhti.com](mailto:info@tamarackhti.com)

ISO13485: 2016 Registered QMS



# GlideWear

## FOOT PROTECTION SOCKS

GlideWear® socks use an **ultra-low friction textile** to protect at-risk areas from the forces that cause diabetic foot ulcers, calluses, blisters, pain & skin breakdown.

GlideWear® is a dual-layer fabric that glides smoothly against itself, **absorbing friction and shear and protecting skin.**

GlideWear® technology is available in four sock models, including forefoot, partial foot, midfoot, and heel/ankle. Each model features SmartKnit® non-binding sock technology.

- ✓ **Ultra-low coefficient of friction**
- ✓ **Clinically effective skin protection**
- ✓ **FDA cleared Class I Medical Device**
- ✓ **Patented, friction-reducing technology**

### SPECIFICATIONS

Sold in Pairs (2/pkg)

Available in all-white or all-black color.

Does not contain latex.

### Adult Shoe Sizes:

	MEN'S	WM'S
<b>S / P</b>	5 - 7	5 - 8
<b>M / M</b>	8 - 10	9 - 11
<b>L / G</b>	11 - 13	12 - 14
<b>XL / XG</b>	14 +	14 +

### CONTENT

75% Polyester  
19% Nylon  
6% Spandex

### CARE INSTRUCTIONS

Wash with mild detergent.



## Forefoot

GW-SX4

### Protection Socks

Indicated primarily for individuals experiencing ulceration or calluses on the toes and metatarsals or forefoot deformities.

Designed for use with orthopedic devices and therapeutic footwear, including walking boots.



## Heel & Ankle

GW-SX6

### Protection Socks

Indicated for individuals at-risk for or experiencing skin breakdown and ulceration of the heel and ankle due to extended bed-stay.

May be used independently or in conjunction with a pressure offloading boot.



## Partial Foot

GW-SX5

### Protection Socks

Indicated for individuals with transmetatarsal (partial foot) amputations to protect the residual foot from breakdown.



## Midfoot

GW-SX8

### Protection Socks

Indicated for individuals with Charcot foot and others at risk for skin breakdown on the midfoot.

Can be used in conjunction with orthotic walking boots (including CROW boots), diabetic shoes, and custom orthotics.

