

Spring Options for MILINE  
Double Action & Dorsiflexion  
Assist Ankle Joints



MILINE Component Size		NO BOOSTER	BOOSTER		
13mm	Stiffness	Low	Moderate	High	—
	Active ROM	15°	25°	12°	—
16mm	Stiffness	Moderate	Moderate	High	Very High
	Active ROM	20°	20°	10°	20°
20mm	Stiffness	Moderate/High	Moderate/High	Very High	Max Stiffness
	Active ROM	15°	15°	10°	15°

\* Spacer used in spring configuration B1