


# FDE ORTHO SPRING CATEGORY (FIRMNESS) SELECTION

## Most Common Starting Point

Category number **3** spring is a good starting point. From here the AFO can be tuned firmer or softer based on the patient needs. The simple four-bolt system makes changing the spring quick and easy.



Categories **.5** and **.75** are not within our common range and are typically used for special situations.

Each patient has unique needs, so there are no specific weight limits or ranges to guide spring selection. Instead, choose the engineered spring category that best enhances the patient's gait and overall functionality.

Categories **6** and **7** are not within our common range and are typically used for special situations.

QUESTIONS? 1.800.FABTECH | FABTECHSYSTEMS.COM