



CD207FK

Fitting Kit

Resilergy™ Dynamic Strut AFO

Weight Range depends on flex
Light to Medium Activity Level



6-month warranty against manufacturer defects, excessive wear or breakage.

External Prosthetic Components



www.coyote.us/instructions



CD207FK.post.revA.11032022




Manufactured by Coyote®
419 N. Curtis Rd., Boise, Idaho 83706 USA
(208) 429-0026 | www.coyote.us



Advena Limited
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta

Fitting Kit Box for the Resilergy™ Dynamic Strut AFO - Map & Parts Included

	More AFO strut flex options will be available			Resilergy Fitting Struts - CD207MS, CD207SS, CD207XS 			AFO Tooling Piece (4) 	AFO Upper Mounting Plate (1) 	AFO Lower Mounting Plate (1) 	AFO Plate Foam (4) 			Button Head Screws 6x25 (2) 	Button Head Screws 6x20 (2) 	Self Leveling Washer Bottom (4) 	Rotation Shim (2) 	Flexion Shim (2) 
---	--	--	---	--	--	---	--	---	---	---	--	---	--	--	--	--	---

EN	Instructions for Use
DE	Gebruiksaanweisung
FR	Notice d'utilisation
ES	Instrucciones para el uso
IT	Istruzioni per l'uso
NO	Bruksanvisning
DA	Bruksanvisning

SV	Bruksanvisning
EL	Οδηγίες Χρήσης
FI	Käyttöohjeet
NL	Gebruiksaanwijzing
PT	Instruções de Utilização
PL	Instrukcja użytkowania
CS	Návod k použití

TR	Kullanım Talimatları
RU	Инструкция по использованию
JA	取扱説明書
ZH	中文说明书
KO	사용 설명서



<https://www.coyote.us/instructions-dynamicstrut>



CD207XS

Extra Soft Flex Kit

Resilergy™ Dynamic Strut AFO

Weight Range 70-150 lbs.
Light to Medium Activity Level



6-month warranty against manufacturer defects, excessive wear or breakage.

External Prosthetic Components



www.coyote.us/instructions



8 11576 02399 3

CD207XS.post.revA.11082022



EC REP

Manufactured by Coyote®
419 N. Curtis Rd., Boise, Idaho 83706 USA
(208) 429-0026 | www.coyote.us



Advena Limited
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta

Parts Included



A. Strut (1) (10.5 x 1.5)

B. Screws (4)

(2) - 6x18 for Top Mounting Plate

(2) - 6x20 for Bottom Mounting Plate

C. Black Screws 6x14 (4)

D. Clamp Washer (4)

E. AFO Mounting Plate (2)

(1) Upper Mounting Plate - Short Posts

(1) Lower Mounting Plate - Tall Posts

F. AFO Fab Tooling Pieces (2)

G. AFO Plate Foam (2)

EN | Instructions for Use

DE | Gebrauchsanweisung

FR | Notice d'utilisation

ES | Instrucciones para el uso

IT | Istruzioni per l'uso

NO | Bruksanvisning

DA | Brugsanvisning

SV | Bruksanvisning

EL | Οδηγίες Χρήσης

FI | Käyttöohjeet

NL | Gebruiksaanwijzing

PT | Instruções de Utilização

PL | Instrukcja użytkowania

CS | Návod k použití

TR | Kullanım Talimatları

RU | Инструкция по использованию

JA | 取扱説明書

ZH | 中文说明书

KO | 사용 설명서



<https://www.coyote.us/instructions-dynamicstrut>



CD207S

Soft Flex Kit

Resilergy™ Dynamic Strut AFO

Weight Range 90-200 lbs.

Light to Medium Activity Level



6-month warranty against manufacturer defects, excessive wear or breakage.

External Prosthetic Components



www.coyote.us/instructions



CD207S.post.revA.11082022



Manufactured by Coyote®
419 N. Curtis Rd., Boise, Idaho 83706 USA
(208) 429-0026 | www.coyote.us



Advena Limited
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta

Parts Included



A. Strut (1) (10.5 x 1.5)

B. Screws (4)

(2) - 6x18 for Top Mounting Plate

(2) - 6x20 for Bottom Mounting Plate

C. Black Screws 6x14 (4)

D. Clamp Washer (4)

E. AFO Mounting Plate (2)

(1) Upper Mounting Plate - Short Posts

(1) Lower Mounting Plate - Tall Posts

F. AFO Fab Tooling Pieces (2)

G. AFO Plate Foam (2)

EN | Instructions for Use

DE | Gebrauchsanweisung

FR | Notice d'utilisation

ES | Instrucciones para el uso

IT | Istruzioni per l'uso

NO | Bruksanvisning

DA | Brugsanvisning

SV | Bruksanvisning

EL | Οδηγίες Χρήσης

FI | Käyttöohjeet

NL | Gebruiksaanwijzing

PT | Instruções de Utilização

PL | Instrukcja użytkowania

CS | Návod k použití

TR | Kullanım Talimatları

RU | Инструкция по использованию

JA | 取扱説明書

ZH | 中文说明书

KO | 사용 설명서



<https://www.coyote.us/instructions-dynamicstrut>



CD207M

Medium Flex Kit

Resilergy™ Dynamic Strut AFO

Weight Range 125-225 lbs.

Light to Medium Activity Level



6-month warranty against manufacturer defects, excessive wear or breakage.

External Prosthetic Components



www.coyote.us/instructions



8 11576 02396 2

CD207M.post.revA.11082022



EC REP

Manufactured by Coyote®
419 N. Curtis Rd., Boise, Idaho 83706 USA
(208) 429-0026 | www.coyote.us



Advena Limited
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta

Parts Included



A. Strut (1) (10.5 x 1.5)

B. Screws (4)

(2) - 6x18 for Top Mounting Plate

(2) - 6x20 for Bottom Mounting Plate

C. Black Screws 6x14 (4)

D. Clamp Washer (4)

E. AFO Mounting Plate (2)

(1) Upper Mounting Plate - Short Posts

(1) Lower Mounting Plate - Tall Posts

F. AFO Fab Tooling Pieces (2)

G. AFO Plate Foam (2)

EN | Instructions for Use

DE | Gebrauchsanweisung

FR | Notice d'utilisation

ES | Instrucciones para el uso

IT | Istruzioni per l'uso

NO | Bruksanvisning

DA | Brugsanvisning

SV | Bruksanvisning

EL | Οδηγίες Χρήσης

FI | Käyttöohjeet

NL | Gebruiksaanwijzing

PT | Instruções de Utilização

PL | Instrukcja użytkowania

CS | Návod k použití

TR | Kullanım Talimatları

RU | Инструкция по использованию

JA | 取扱説明書

ZH | 中文说明书

KO | 사용 설명서



<https://www.coyote.us/instructions-dynamicstrut>