


## Secure suspension made easy.

Safe, predictable alignment and suspension. **Just turn the dial.**

### How it works

#### Don:

Snap into your liner and turn the dial to gently pull your limb into your socket. (Goodbye skin shear.)




#### Adjust:

Adjust your suspension by turning the dial. (No need to don and doff.)



#### Doff:

Turn the dial to release your liner, then press the quick-connect snap to disconnect.



*You can don and doff standing or sitting.*

"The RevoLock<sup>®</sup> system lets me dial myself down into my socket, so I have the **right fit** and suspension **every time**."

- Kim "Tin Man" Acheson, RevoLock user



### Who is it for?

If any of these sound like you, **RevoLock<sup>®</sup>** could be your perfect solution.

#### Do you:

- Struggle with your current suspension system?
- Have trouble locating a pin lock?
- Have long or bulbous limbs?
- Have large volume changes during the day?
- Struggle with balance when donning your socket?
- Have a limited range of motion?

**Talk to your practitioner about how RevoLock<sup>®</sup> can help you focus more on your life and less on your socket.**

To find out how adjustability can work for you, connect with a user and more, scan the QR code.



**Let's  
Improve  
Lives**

Click Medical's mission is to empower users to overcome limitations, regain mobility, and reclaim independence.

[www.clickmedical.co](http://www.clickmedical.co)