



All-Purpose Boot

All season closed toe protection. Deep opening allows for bulky dressing. Heel locking ankle strap and breathable padded upper. May also be used with synthetic casts. Sold each.



Open Toe Cast Boot

Economical canvas cast boot with three hook and loop closures. Triple layer EVA foam rocker bottom provides shock absorption and eases roll into and out of the gait cycle. Sold each.



Med-Surg™ Shoe

Square toe design acts as a bumper to provide additional protection. Rocker sole reduces plantar pressure. Also available in Original style shown in photo above. Sold each.



OrthoWedge™ Shoe

Lessens friction and pressure by reducing weight bearing under the forefoot. Ideal following surgery or as a healing shoe for diabetic or abnormal foot conditions. Sold each.



Softie™ Shoe

Economical, post-op shoe for soft tissue trauma or surgery. Semi-rigid metatarsal shank for forefoot control. Conformable, breathable upper material. Easy closure. Sold each.



Wound Care Shoe

The ultimate shoe for healing diabetic foot ulcers and infections. Deep rocker sole contains 4 layers of differing density insoles. Uppers have removable sections. Sold in pairs.



HeelWedge™ Shoe

Protects the heel by reducing weight bearing pressure. Locking ankle strap to assure secure and safe ambulation. Used to off-load ulcers, infections or trauma to the heel. Sold each.



SlimLine™ Cast Boot

Designed specifically for the thinner synthetic cast. Incorporates dual rocker sole with custom fitted appearance. Stretchable upper keeps foot seated in insole. Sold each.



DARCO® HEALING SHOES — Continued

All-Purpose Boot	X-Small	Small	Medium	Large	X-Large
Men's	3-4.5	5-6.5	7-8.5	9-11	11.5-14
Women's	3-6.5	7-9.5	10-12.5	13+	-----
Black	APQ0B	APQ1B	APQ2B	APQ3B	APQ4B

Open Toe Cast Boot	Small	Medium	Large	X-Large
Men's	-----	6-8.5	9-11.5	12-14
Women's	5.5-7.5	8-10	-----	-----
Blue	CBO1	CBO2	CBO3	CBO4

Med-Surg™ Shoe	Small	Medium	Large	X-Large
Men's	6-8	8.5-10	10.5-12	12.5-14
Black	MQM1B	MQM2B	MQM3B	MQM4B
Navy (Original)	MSM1N	MSM2N	MSM3N	MSM4N
Women's	4-6	6.5-8	8.5-10	-----
Black	MQW1B	MQW2B	MQW3B	
Navy (Original)	MSW1N	MSW2N	MSW3N	

OrthoWedge™ Shoe	X-Small	Small	Medium	Large	X-Large
Men's	-----	6-8	8.5-10	10.5-12	12.5-14
Women's	4-7	7.5-10	-----	-----	-----
Black	OQ0B	OQ1B	OQ2B	OQ3B	OQ4B

Softie™ Shoe	Small	Medium	Large	X-Large
Men's	6-8	8.5-10	10.5-12	12.5-14
Black	STM1B	STM2B	STM3B	STM4B
Women's	4-6	6.5-8	8.5-10	-----
Black	STW1B	STW2B	STW3B	

Wound Care Shoe System	Small	Sm/Med	Medium	Large	X-Large	XX-Large
Men's	6-7.5	8-8.5	9-9.5	10-10.5	12-13.5	14-16
Women's	6-8.5	9-9.5	10-10.5	11-12.5	13-14.5	15-17
Black	WCS1B	WCS7B	WCS2B	WCS3B	WCS4B	WCS5B

HeelWedge™ Shoe	X-Small	Small	Medium	Large	X-Large
Men's	-----	6-8	8.5-10	10.5-12	12.5-14
Women's	4-7	7.5-10	-----	-----	-----
Black	HQ0B	HQ1B	HQ2B	HQ3B	HQ4B

SlimLine™ Cast Boot	X-Small	Small	Medium	Large	X-Large
Men's	-----	-----	6-8	8.5-11	11.5-14
Women's	3-5	5.5-8	8.5-10	-----	-----
Navy	SLO0N	SLO1N	SLO2N	SLO3N	SLO4N