Symes Prosthetic Sock


STANDARD SYMES SOCK SIZES — IN SUPER SOCK AND SOFT-SOCK

The Standard Size Symes Socks listed below will provide a well-fitting sock for most amputees—if special sizes are needed, please refer to pages 34 and 35. Sold each.

Appro	Approximate Range—					
Size	Ankle Width	Ball Width	Ankle Circumference	Ball Circumference		
А	3" (7.5cm)	3.50" (9cm)	5"-6.50" (12.5-16.5cm)	6.50"-8" (16.5-20cm)		
В	3" (7.5cm)	4" (10cm)	5"-6.50" (12.5-16.5cm)	7.50"-9.50" (19-22.5cm)		
С	3.50" (9cm)	4" (10cm)	6"-8" (15-20cm)	8.50"-9.50" (20-24cm)		
D	3.50" (9cm)	4.50" (11cm)	6"-8" (15-20cm)	8.50"-10" (21-25cm)		
Е	4" (10cm)	4.50" (11cm)	7.50"-10.50" (19-26.5cm)	9-10.50" (22.5-26.5cm)		
Н	5" (12.5cm)	5.50" (14cm)	9.50"-11" (24-28cm)	10.50"-12" (26.5-30.5cm)		

The toe width at the very tip end of Symes Socks is wider (on both standard Symes sizes and special Symes sizes) than the width at the very tip end of regular B/K size prosthetic socks. Use of regular B/K size socks for Symes amputees may result in wrinkles at the bottom of the limb and/or break through the sock near the toe.

STANDARD SYMES SOCK — SUPER SOCK



Example of sizing 1S15-AO-18: 18" x 6" x 3" x 3.50" Graduated for larger sizes.

TOP X ANKLE 6" x 3" x 3.50 (15 x 7.5 x 90 3-PLY)"	6" x 3" x 4" (15 x 7.5 x 10 3-PLY	OCM) 5-PLY	6" x 3.50" x 4 (15 x 9 x 10C 3-PLY		6" x 3.50" x 4 (15 x 9 x 11C 3-PLY		LENGTH
Super Sock —	- Size O							
1S13-A0-18	1S15-A0-18	1S13-B0-18	1S15-B0-18	1S13-C0-18	1S15-C0-18	1S13-D0-18	1S15-D0-18	18" (45cm)
1S13-A0-20	1S15-A0-20	1S13-B0-20	1S15-B0-20	1S13-C0-20	1S15-C0-20	1S13-D0-20	1S15-D0-20	20" (50cm)
1S13-A0-22	1S15-A0-22	1S13-B0-22	1S15-B0-22	1S13-C0-22	1S15-C0-22	1S13-D0-22	1S15-D0-22	22" (55cm)
1S13-A0-24	1S15-A0-24	1S13-B0-24	1S15-B0-24	1S13-C0-24	1S15-C0-24	1S13-D0-24	1S15-D0-24	24" (60cm)
TOP X ANKLE X BALL 7" x 3" x 3.50" (17.5 x 7.5 x 9CM)		7" x 3" x 4" (17.5 x 7.5 x	10CM)	7" x 3.50" x 4 (17.5 x 9 x 10		7" x 3.50" x 4 (17.5 x 9 x 11		
Super Sock —	- Size 1							
1513-A1-18	1S15-A1-18	1S13-B1-18	1S15-B1-18	1S13-C1-18	1S15-C1-18	1S13-D1-18	1S15-D1-18	18" (45cm)
1513-A1-20	1S15-A1-20	1S13-B1-20	1S15-B1-20	1513-C1-20	1S15-C1-20	1S13-D1-20	1S15-D1-20	20" (50cm)
1S13-A1-22	1S15-A1-22	1S13-B1-22	1S15-B1-22	1S13-C1-22	1S15-C1-22	1513-D1-22	1S15-D1-22	22" (55cm)
1S13-A1-24	1S15-A1-24	1S13-B1-24	1S15-B1-24	1S13-C1-24	1S15-C1-24	1\$13-D1-24	1\$15-D1-24	24" (60cm)

STANDARD SYMES SOCK — SOFT-SOCK

TOP x ANKLE 6" x 3" x 3.50 (15 x 7.5 x 90)"			7" x 4" x 4.50 (18 x 10 x 110		7" x 5" x 5.50 (18 x 12.5 x 1		
Super Sock —	- Size 0			Soft-Sock — S	ize 1			
1SP3-A0-18	1SP5-A0-18	_	—	1SP3-E1-18	1SP5-E1-18	1SP3-H1-18	1SP5-H1-18	18" (45cm)
1SP3-A0-20	1SP5-A0-20	_	_	1SP3-E1-20	1SP5-E1-20	1SP3-H1-20	1SP5-H1-20	20" (50cm)
1SP3-A0-22	1SP5-A0-22	_	_	1SP3-E1-22	1SP5-E1-22	1SP3-H1-22	1SP5-H1-22	22" (55cm)
1SP3-A0-24	1SP5-A0-24	—	—	1SP3-E1-24	1SP5-E1-24	1SP3-H1-24	1SP5-H1-24	24" (60cm)

Prosthetic Sock Ordering Guidelines

KNIT-RITE STANDARD PROSTHETIC SOCK SIZING

Standard Sizing Guidelines*		
Knit-Rite Prosthetic Sock	Top Width	Toe Width
Size No. 0	6" (15cm)	3.50" (9cm)
Size No. 1	7" (17.5cm)	4" (10cm)
Size No. 2	8" (20cm)	5" (12.5cm) [over 24" (60cm) length 4" (10cm)]
Size No. 3	9" (23cm)	6" (15cm) [over 20" (50cm) length 5.50" (14cm)]
Size No. 4	10" (25cm)	7" (17.5cm) [over 20" (50cm) length 6" (15cm)]
Size No. 5	11" (27.5cm)	7" (17.5cm) [over 18" (45cm) length 5.50" (14cm)]
Size No. 6	12" (30cm)	8" (20cm) [over 18" (45cm) length 6" (15cm)]

Knit-Rite Arm Sock	Top Width	Toe Width
Knit-Rite Size A; VA Size Arm 0	4" (10cm)	3" (7.5cm)
Knit-Rite Size B; VA Size Arm 1	5.50" (14cm)	3.50" (9cm)

*These sizes shipped unless otherwise specified. Sterling Prosthetic Socks have some variation from Knit-Rite sizing.

Also see pages for Standard PTB Sizes, Symes, Hip Socks, X-Wool, Soft-Sock and A-Plus.

HCPCS Bi	illing Codes
Codes	Descriptions
L0984	Protective Body Sock, each
L2820	Addition to lower extremity orthosis soft interface for molded plastic, below the knee section (Knit-Rite AFO Liner Sock).
L2840	Addition to lower extremity orthosis, tibial length sock, fracture or equal, each
L2850	Addition to lower extremity orthosis, sock, femoral length sock, fracture or equal, each
L3995	Addition to upper extremity orthosis, sock, fracture or equal, each
L8400	Prosthetic Sheath, below knee, each
L8410	Prosthetic Sheath, above knee, each
L8415	Prosthetic Sheath, upper limb, each
L8420	Prosthetic Sock, multiple ply, below knee, each*
L8430	Prosthetic Sock, multiple ply, above knee, each*
L8435	Prosthetic Sock, multiple ply, upper limb, each* *Traditional Knit-Rite Prosthetic Socks with knit-in taper, including Knit-Rite Liner Liner. Multiple (1-8) plies of natural or synthetic yarns. Knit on flat bed knitting machines.
L8440	Prosthetic Shrinker, below knee, each
L8460	Prosthetic Shrinker, above knee, each
L8465	Prosthetic Shrinker, upper limb, each
L8470	Fitting Sock, single ply, below knee, each
L8480	Fitting Sock, single ply, above knee, each
L8485	Fitting Sock, single ply, upper limb, each Fitting Socks, single ply, tubular construction. Knit on circular knitting machines.
L7600	Prosthetic Donning Sleeve, any material, each

We have included the above codes to assist the provider in determining and obtaining proper reimbursement for our products. It is the provider's responsibility to ascertain the appropriateness of individual applications of these billing codes.