

# Townsend XLR8

Pneumatic Walker

Available In 5 Sizes, Full  
Height Or Mid Calf



# XLR8 Technology



## Anatomically Correct Calf Section

The upper edge of the posterior calf shell is flared to enhance comfort and patient compliance.



## Individual Medial & Lateral Built-In Air Pumps

User-friendly air pumps enable individual adjustment of compression air chambers on both sides of the ankle.



## Adjustable Front Plate Compression (Full Height Model)

The lower leg strap can be inserted through slots in both sides of the posterior shell to create additional anterior and circumferential compression.

## Three Toe Guard Settings

Based on the patient's actual foot length, the adjustable toe guard can be repositioned or removed.

# Townsend XLR8 Walker

The Townsend XLR8 (Accelerate) is a full shell pneumatic walking boot designed to enable patients to safely and comfortably ambulate after an injury or surgery. The XLR8 Walker features integrated medial and lateral air pumps (and air release valves) for individual control of the dual air chambers. The built-in blue pumps are easily accessible and a few quick squeezes will inflate the heavy duty air chambers that are concealed between the rigid outer shell and soft foam liner. This anterior-closure liner has adjustable sizing and extends two inches above the upper shell configuration to eliminate potential edge pressure or soft tissue impingement.

The XLR8 has a soft internal foot pad and wide heel area to accommodate bandages, swelling and inflation of the air chambers. The XLR8 also features a non-skid, low profile, shock-absorbing outer sole with a rocker motion that allows for smooth transition from heel strike to midstance and toe off. The wide base and tread pattern is designed for durability as well as variety of surfaces and weather conditions.



Ordered by: \_\_\_\_\_ Phone # ( \_\_\_\_\_ ) \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_ Townsend Account # \_\_\_\_\_

Bill To: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: ( \_\_\_\_\_ ) \_\_\_\_\_ Fax: ( \_\_\_\_\_ ) \_\_\_\_\_

Phone: ( \_\_\_\_\_ ) \_\_\_\_\_ Fax: ( \_\_\_\_\_ ) \_\_\_\_\_

Shipping Preference: ☐ Ground ☐ 2-Day P.M. ☐ 2-Day A.M. ☐ Next Day P.M. ☐ Next Day A.M.  
(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

## Full Height Product Codes:

- ☐ **24-113 XSmall** -- Quantity: \_\_\_\_\_
- ☐ **24-104 Small** -- Quantity: \_\_\_\_\_
- ☐ **24-105 Medium** -- Quantity: \_\_\_\_\_
- ☐ **24-106 Large** -- Quantity: \_\_\_\_\_
- ☐ **24-114 XLarge** -- Quantity: \_\_\_\_\_

## Mid Calf Product Codes:

- ☐ **24-115 XSmall** -- Quantity: \_\_\_\_\_
- ☐ **24-116 Small** -- Quantity: \_\_\_\_\_
- ☐ **24-117 Medium** -- Quantity: \_\_\_\_\_
- ☐ **24-118 Large** -- Quantity: \_\_\_\_\_
- ☐ **24-119 XLarge** -- Quantity: \_\_\_\_\_

## Sizing Reference

Item #	Size	Men	Women
24-113	XSmall	Up to 4	Up to 5
24-104	Small	4 - 7	5 - 8
24-105	Medium	7 - 10	8 - 11
24-106	Large	10 - 13	11 - 15
24-114	XLarge	13 +	15 +

You can photocopy this page and use it to fax or email XLR8 orders to Townsend. The fax number is 800.798.2722 and you can email scanned order forms to [eorders@townsenddesign.com](mailto:eorders@townsenddesign.com). If you prefer to call in your order, the direct line to Customer Service is 800.700.2722, and customer representatives are routinely available to provide assistance between 6:00 a.m. and 4:30 p.m. The information on this page will enable you to reference the product codes and other information that will need to be provided during phone orders.

