

Extension Adapters

When additional height is needed in a prosthesis, extension adapters can be valuable components. Available in various heights, extension adapters are offered with pyramids and receivers.

HEIGHT EXTENSION ADAPTERS



Great for use with Fusion or DuraLite™ Feet.

Height Extension Adapters are ideal for when additional height is needed in a prosthesis. Applications may include adding height when a shank-style foot is trimmed or when other components are no longer desired in a prosthesis. Available in 24 mm and 46 mm. The titanium adapters have a pyramid on one end and pyramid receiver on the other.

PRODUCT WEIGHT 46 mm: 96 g
24 mm: 79 g

PATIENT WEIGHT LIMIT* 300 lbs (136 kg)

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
46 mm Height Extension Adapter	FND-136146
24 mm Height Extension Adapter	FND-136124

TITANIUM DUAL PYRAMID RECEIVERS



With pyramid receivers on both ends, these titanium adapters provide excellent adjustability options when clearance is limited. They are available in 32 mm, 45 mm, 60 mm, and 75 mm heights.

PRODUCT WEIGHT 75 mm: 126 g 45 mm: 92 g
60 mm: 107 g 32 mm: 82 g

PATIENT WEIGHT LIMIT* 350 lb (160 kg) for U.S. Levels 1-3; 300 lb (136 kg) for U.S. Level 4

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
75 mm Titanium Dual Pyramid Receiver	FND-228075
60 mm Titanium Dual Pyramid Receiver	FND-228060
45 mm Titanium Dual Pyramid Receiver	FND-228045
32 mm Titanium Dual Pyramid Receiver	FND-228032

MAGNUM TITANIUM DOUBLE ENDED PYRAMID



For use with adult-sized pyramid receivers, this component features pyramids on each end. May be used in heavy duty applications.

PRODUCT WEIGHT 66 g

PATIENT WEIGHT LIMIT* 350 lb (160 kg)

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
Magnum Titanium Double Ended Pyramid	MHD-DE-P

*Body weight plus any loads normally or routinely carried cannot exceed this weight limit.