

# Titanium 4-Hole Pyramid Adapters

## TITANIUM ROTATABLE PYRAMID ADAPTER



The Titanium Rotatable Pyramid Adapter is lightweight and provides 360 degrees of rotation. Set screws are included.

**PRODUCT WEIGHT** 57 g

**PATIENT WEIGHT LIMIT\*** 250 lb (115 kg)

**WARRANTY** 12 months from date of patient fitting

### ITEM NUMBERS

Product	Item #
Titanium Rotatable Pyramid Adapter	FND-228003

## TITANIUM 4-HOLE PYRAMID ADAPTER



The 4-Hole Titanium Pyramid Adapter attaches to any socket or knee with a distal 4-hole pattern.

**PRODUCT WEIGHT** 56 g

**PATIENT WEIGHT LIMIT\*** 250 lb (115 kg)

**WARRANTY** 12 months from date of patient fitting

### ITEM NUMBERS

Product	Item #
Titanium 4-Hole Pyramid Adapter	FND-227014

## TITANIUM 4-HOLE PYRAMID, 20 MM OFFSET



The Titanium 4-Hole Pyramid Adapter, 20 mm Offset attaches to any socket or knee with a distal 4-hole pattern. The pyramid is offset 20 mm anterior/posterior or medial/lateral from the center. Screws are included.

**PRODUCT WEIGHT** 93 g

**PATIENT WEIGHT LIMIT\*** 250 lb (115 kg)

**WARRANTY** 12 months from date of patient fitting

### ITEM NUMBERS

Product	Item #
Titanium 4-Hole Pyramid, 20 mm Offset	FND-224120

\*Body weight plus any loads normally or routinely carried cannot exceed this weight limit.

**TITANIUM 4-HOLE SLIDE ADAPTER WITH ROTATION**

The Titanium 4-Hole Slide Adapter with Rotation provides 360 degrees of rotation and 1/2" of slide. Set screws are included.

**PRODUCT WEIGHT** 125 g

**PATIENT WEIGHT LIMIT\*** 350 lb (160 kg) for U.S. Levels 1-3; 300 lb (136 kg) for U.S. Level 4

**WARRANTY** 12 months from date of patient fitting

**ITEM NUMBERS**

Product	Item #
Titanium 4-Hole Slide Adapter with Rotation	FND-224150

**LIMBLOGIC® VACUUM PYRAMID**

The LimbLogic Vacuum Pyramid is suitable for instances when there is insufficient clearance to mount the LimbLogic pump directly below the socket. The Vacuum Pyramid allows practitioners to relocate the controller while maintaining the in-line design of a prosthesis. Additionally, practitioners may use the Vacuum Pyramid to retrofit existing vacuum sockets to use with LimbLogic.

**PRODUCT WEIGHT** 62 g

**PATIENT WEIGHT LIMIT\*** 350 lb (160 kg) for U.S. Levels 1-3; 300 lb (136 kg) for U.S. Level 4

**WARRANTY** 12 months from date of patient fitting

**ITEM NUMBERS**

Product	Item #
LimbLogic® Vacuum Pyramid	LLV-01044