



The ALPS™ Skin Reliever eliminates shear forces, friction, and prevents abrasions on the skin of the stump. It has both superior comfort and durability combined. It also accommodates volume fluctuations over the life of the prosthesis. There is superior elongation for very little compression against the skin. There is no need to make a new socket, ALPS Skin Reliever can be worn under silicone, gel, and pelite liners.

- Replaces pelite and socks
- High Elasticity
- Reduces shearing and abrasions
- Softens and soothes skin

PRODUCT NO.

ENCP-XX


Skin Reliever, Tapered from 3mm - 1.5mm

XX indicates size, see sizing chart and illustrations on reverse for details on sizing and differences between configurations.

Patent Pending



Skin Reliever ENCP Cross Section Diagrams
3mm - 1.5mm Tapered - No Fabric or Distal Attachment



Skin Reliever ENCP Sizing Chart
Measurements taken 6cm from the distal end

Measured Size (cm) ^(*)	ENCP
Description	Tapered
10-19 cm	ENCP-16
20-27 cm	ENCP-20
28-31 cm	ENCP-28

How to Size
Measurements taken 6cm from the distal end

