



The Padsil was developed to increase patient comfort in combination with ALPS™ liners. By increasing distal pressure and enhancing vascular function in distal tissues, the more active user will appreciate the additional cushioning, softness, and durability. Users who experience volume changes will also appreciate the Padsils ability to lift the limb back into the proper position. They can be placed at the bottom of the socket or inside a pin locking liner.

- Soft medical grade silicone
- Keeps its shape over time
- Helps maintain total contact
- Can be modified
- Blue dummy padsils are available for laminating

### PRODUCT NO.

PDS-XX


**Silicone Distal Pads**

XX indicates size, see sizing chart and illustrations on reverse for details on sizing and differences between configurations.



Padsil Fitting Diagrams

Can be used in any system



Padsil Sizing Chart

Measurements taken 6cm from the distal end

Measured Size (cm) <sup>(*)</sup>	Padsil
Description	Size
14-19 cm	PDS-S
20-23 cm	PDS-M
24-27 cm	PDS-L
28-31 cm	PDS-XL
32-35 cm	PDS-XXL

How to Size

Measurements taken 6cm from the distal end

