

Prosthetic Lubricant

PA604



70% Isopropyl Alcohol Prosthetic Lubricant for Liners

ALPS™ Prosthetic Lubricant is used to help in donning prosthetic liners. It allows the patient to don the liner or sheath easily while eliminating the sticky feel.

Recommended for use with:

- ☒ Gel (Thermoplastic) - All ALPS™ Gel Liners
- ☒ Silicone - CLEARPRO
- ☒ Urethane

PRODUCT NUMBER

PA604 **Prosthetic Lubricant**, 15 oz.

PA604 BX12 **Prosthetic Lubricant**, 15 oz. Box of 12



CAN#6863 Form #100644 Rev B



Tel: 727.528.8566 | Tel: 800.574.5426
www.easyliner.com | info@easyliner.com