



**NEW  
LAMINATION  
ADAPTERS  
AND  
FLEXION  
PLATE  
SEE PAGE 6 & 8**

# ENDOSKELETAL

*Endoskeletal Components*




**2019**

CATALOG



## NEW LAMINATION ADAPTERS


College Park now offers a robust line of lamination adapters! In our efforts to diversify our line of high-quality endo components, we now carry seven new lamination plates as well as one new flexion plate. This new selection accommodates almost any socket design and is made from materials such as aluminum, aircraft-grade titanium, and heavy-duty stainless steel.



## NOISE ABATEMENT

90% of noise occurs at the interfaces, and there are 19 mating surfaces in a typical prosthetic leg. College Park utilizes premium materials and special treatments for the screws and threads that can make all the difference. Designed and manufactured by College Park in the USA, the benefits of premier endo components are available at no additional cost to you.

When two mating surfaces are smooth, the potential for noise due to friction is almost guaranteed. Contact treatment and integrated ridges on our pyramids assist in reducing this micro-movement for noise abatement right at the clamping surface.



## SUPER SET SCREWS

Our steel set screws are manufactured for unique benefits not found elsewhere. They are roughened to provide noise abatement against the mating surface, then treated with a zinc coating. This corrosion-resistant coating protects against harsh elements to increase the lifespan of the screws.

The set screws are also treated with a nylon patch. The nylon brings a serviceable, reusable alternative to Loctite. There's no more forgetting to apply or over-applying. It prevents the screws from backing out, yet won't strip the threads when removing.



518 LIT CAT END 190426

©2019, College Park Industries, Inc. All rights reserved.



# CONTENTS

## 2019

### COMPONENTS

<i>Titanium</i>	<i>pg. 2</i>
<i>Heavy Duty</i>	<i>pg. 3</i>
<i>Economy</i>	<i>pg. 3</i>
<i>Pediatric</i>	<i>pg. 4</i>
<i>Geriatric</i>	<i>pg. 4</i>
<i>Kits</i>	<i>pg. 5</i>

### LAMINATION

<i>Adapters</i>	<i>pg. 6</i>
<i>Kits</i>	<i>pg. 7</i>
<i>Additional Parts</i>	<i>pg. 7</i>

### FLEXION PLATES

*pg. 8*

### ADDITIONAL HARDWARE

<i>Standard</i>	<i>pg. 8</i>
<i>Pediatric</i>	<i>pg. 8</i>



**college park**

TECHNOLOGY *for the* HUMAN RACE


# TITANIUM

components

**MATERIAL**  
TITANIUM

**WEIGHT LIMIT**  
300 lbs (136 kg)

COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>LOW-PROFILE TUBE ADAPTER</b> 250 mm 30 mm diameter	57 mm (min) 221 g (max)	<b>TA T250 30</b>
 <b>LOW-PROFILE TUBE ADAPTER</b> 430 mm 30 mm diameter	57 mm (min) 337 g (max)	
 <b>TUBE CLAMP</b> 30 mm diameter	51 mm 82 g	<b>TC T 30</b>
 <b>DOUBLE ENDED ADAPTER,</b> 35 mm	35 mm 95 g	
 <b>DOUBLE ENDED ADAPTER</b> 45 mm	45 mm 110 g	<b>DR T45</b>
 <b>DOUBLE ENDED ADAPTER</b> 60 mm	60 mm 123 g	
 <b>DOUBLE ENDED ADAPTER</b> 75 mm	75 mm 136 g	<b>DR T75</b>




COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>4-HOLE PYRAMID RECEIVER</b> Standard	17 mm 64 g	<b>4HPR T</b>
 <b>4-HOLE PYRAMID RECEIVER</b> Rotatable	19 mm 90 g	
 <b>4-HOLE PYRAMID ADAPTER</b> Standard	7 mm 6 g	<b>4HPA T</b>
 <b>4-HOLE PYRAMID ADAPTER</b> Rotatable	9.5 mm 77 g	
 <b>4-HOLE PYRAMID OFFSET</b> Axial	7 mm 55 g	<b>4HPA AO T</b>
 <b>4-HOLE PYRAMID OFFSET</b> Diagonal	7 mm 64 g	
 <b>4-HOLE PYRAMID ADAPTER</b> w/Pin-Lock Hole	7 mm 44 g	<b>4HPH T</b>
 <b>4-HOLE PYRAMID ADAPTER</b> w/Pin-Lock Hole Rotatable	7 mm 74 g	



# HEAVY DUTY

## components

**MATERIAL**  
STAINLESS STEEL

**WEIGHT LIMIT**  
400 lbs (181 kg)

COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>LOW-PROFILE TUBE ADAPTER</b> 250 mm 34 mm diameter	57 mm (min)	<b>TA S250 34</b>
	302 g (max)	
 <b>LOW-PROFILE TUBE ADAPTER</b> 430 mm 34 mm diameter	57 mm (min)	<b>TA S430 34</b>
	425 g (max)	
 <b>TUBE CLAMP</b> 34 mm diameter	55 mm	<b>TC S 34</b>
	146 g	




COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>4-HOLE PYRAMID ADAPTER</b> Standard	7 mm	<b>4HPA S</b>
	118 g	
 <b>4-HOLE PYRAMID RECEIVER</b> Standard	17 mm	<b>4HPR S</b>
	75 g	



# ECONOMY

## components

**MATERIAL**  
STAINLESS STEEL

**WEIGHT LIMIT**  
300 lbs (136 kg)

COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>LOW-PROFILE TUBE ADAPTER</b> 250 mm 30 mm diameter	57 mm (min)	<b>TA S250 30</b>
	303 g (max)	
 <b>LOW-PROFILE TUBE ADAPTER</b> 430 mm 30 mm diameter	57 mm (min)	<b>TA S430 30</b>
	418 g (max)	
 <b>TUBE CLAMP</b> 30 mm diameter	51 mm	<b>TC S 30</b>
	140 g	

COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>4-HOLE PYRAMID ADAPTER*</b> Standard	7 mm	<b>4HPA S</b>
	118 g	
 <b>4-HOLE PYRAMID RECEIVER*</b> Standard	17 mm	<b>4HPR S</b>
	75 g	

# PEDIATRIC

components

**MATERIAL**  
TITANIUM

**WEIGHT LIMIT**  
up to 165 lbs (75 kg)

## ADDITIONAL HARDWARE

COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>LOW-PROFILE TUBE ADAPTER*</b> 300 mm 22 mm diameter	47 mm (min) 169 g (max)	<b>TA T300 22</b>
 <b>TUBE CLAMP*</b> 22 mm diameter	39 mm 46 g	<b>TC T 22</b>
 <b>LOW-PROFILE TUBE ADAPTER</b> 300 mm 25 mm diameter	47 mm (min) 209 g (max)	<b>TA T300 25</b>
 <b>TUBE CLAMP</b> 25 mm diameter	39 mm 45 g	<b>TC T 25</b>
 <b>4-HOLE PYRAMID RECEIVER</b> Standard	15 mm 44 g	<b>4HPR TP</b>
 <b>4-HOLE PYRAMID RECEIVER</b> Rotatable	16 mm 50 g	<b>4HPRR TP</b>
 <b>4-HOLE PYRAMID ADAPTER</b> Standard	8 mm 36 g	<b>4HPA TP</b>
 <b>4-HOLE PYRAMID ADAPTER</b> Rotatable	9 mm 56 g	<b>4HPAR TP</b>



HARDWARE	CONTENTS	PART NUMBER
 <b>4-HOLE FASTENER KIT</b> (4pc) M5 x 14 mm, Pediatric	Four fasteners (14 mm)	<b>4HFK 14P</b>
 <b>SUPER SET SCREWS</b> (8pc) M6 x 12 mm, Pediatric	Eight Super Set Screws	<b>CP SSET P</b>
 <b>TUBE CLAMP SCREW</b> M4 x 14 mm, Pediatric	One tube clamp screw	<b>CP TCS P</b>

# GERIATRIC

components

**MATERIAL**  
TITANIUM

**WEIGHT LIMIT**  
165 lbs (75 kg)

COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>LOW-PROFILE TUBE ADAPTER</b> 300 mm 25 mm diameter with 30 mm adult receiver	52 mm (min) 242 g (max)	<b>TA T300 2530</b>
 <b>TUBE CLAMP</b> 25 mm diameter with 30 mm adult receiver	45 mm 83 g	<b>TC T 2530</b>

# COMPONENT kits

TITANIUM  
ECONOMY  
HEAVY DUTY  
PEDIATRIC  
GERIATRIC

## TITANIUM

KIT	INCLUDED COMPONENTS	PART NUMBER
<b>TRANSTIBIAL KIT W/ PYRAMID ADAPTER</b> Titanium	<ul style="list-style-type: none"> <li>• 30 mm Tube Adapter, 250 mm</li> <li>• 30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter</li> </ul>	<b>TI TTK PA</b>
<b>TRANSTIBIAL KIT W/ PYRAMID ADAPTER WITH PIN-LOCK HOLE</b> Titanium	<ul style="list-style-type: none"> <li>• 30 mm Tube Adapter, 250 mm</li> <li>• 30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter (w/pin-lock hole)</li> </ul>	<b>TI TTK PAH</b>

KIT	INCLUDED COMPONENTS	PART NUMBER
<b>TRANSFEMORAL KIT W/ PYRAMID ADAPTER</b> Titanium	<ul style="list-style-type: none"> <li>• 30 mm Tube Adapter, 430 mm</li> <li>• 30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter</li> </ul>	<b>TI TFK PA</b>
<b>TRANSFEMORAL KIT W/ PYRAMID ADAPTER WITH PIN-LOCK HOLE</b> Titanium	<ul style="list-style-type: none"> <li>• 30 mm Tube Adapter, 430 mm</li> <li>• 30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter (w/pin-lock hole)</li> </ul>	<b>TI TFK PAH</b>

## ECONOMY

KIT	INCLUDED COMPONENTS	PART NUMBER
<b>TRANSTIBIAL KIT W/ PYRAMID ADAPTER</b> Stainless Steel	<ul style="list-style-type: none"> <li>• 30 mm Tube Adapter, 250 mm</li> <li>• 30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter</li> </ul>	<b>SS TTK PA</b>
<b>TRANSFEMORAL KIT W/ PYRAMID ADAPTER</b> Stainless Steel	<ul style="list-style-type: none"> <li>• 30 mm Tube Adapter, 430 mm</li> <li>• 30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter</li> </ul>	<b>SS TFK PA</b>

## HEAVY DUTY

KIT	INCLUDED COMPONENTS	PART NUMBER
<b>TRANSTIBIAL KIT W/ PYRAMID ADAPTER</b> Stainless Steel	<ul style="list-style-type: none"> <li>• 34 mm Tube Adapter, 250 mm</li> <li>• 34 mm Tube Clamp</li> <li>• 4-hole pyramid adapter</li> </ul>	<b>HD TTK PA</b>
<b>TRANSFEMORAL KIT W/ PYRAMID ADAPTER</b> Stainless Steel	<ul style="list-style-type: none"> <li>• 34 mm Tube Adapter, 430 mm</li> <li>• 34 mm Tube Clamp</li> <li>• 4-hole pyramid adapter</li> </ul>	<b>HD TFK PA</b>

## PEDIATRIC & GERIATRIC

KIT	INCLUDED COMPONENTS	PART NUMBER
<b>22mm KIT W/ PEDIATRIC PYRAMID ADAPTER</b> Titanium	<ul style="list-style-type: none"> <li>• 22 mm Tube Adapter, 300 mm</li> <li>• 22 mm Tube Clamp</li> <li>• 4-hole pyramid adapter (pediatric)</li> </ul>	<b>22K PAP</b>
<b>25mm KIT W/ PEDIATRIC PYRAMID ADAPTER</b> Titanium	<ul style="list-style-type: none"> <li>• 25 mm Tube Adapter, 300 mm</li> <li>• 25 mm Tube Clamp</li> <li>• 4-hole pyramid adapter (pediatric)</li> </ul>	<b>25K PAP</b>


KIT	INCLUDED COMPONENTS	PART NUMBER
<b>25-30mm KIT W/ PEDIATRIC PYRAMID ADAPTER</b> Titanium	<ul style="list-style-type: none"> <li>• 25-30 mm Tube Adapter, 300 mm</li> <li>• 25-30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter (pediatric)</li> </ul>	<b>2530K PAP</b>
<b>25-30mm KIT W/ ADULT PYRAMID ADAPTER</b> Titanium	<ul style="list-style-type: none"> <li>• 25-30 mm Tube Adapter, 300 mm</li> <li>• 25-30 mm Tube Clamp</li> <li>• 4-hole pyramid adapter (adult)</li> </ul>	<b>2530K PA</b>


# LAMINATION

*adapters*

**MATERIAL**  
MULTI

**WEIGHT LIMIT**  
300 lbs (136 kg)

COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>4-EAR SOCKET ADAPTER W/ PIN LOCK HOLE</b> Titanium	5.7 mm 55.2 g	<b>4ESAH T</b>
 <b>4-EAR SOCKET ADAPTER, W/ PIN LOCK HOLE</b> Stainless Steel	5.7 mm 96.5 g	<b>4ESAH S</b>
 <b>THREADED ADAPTER W/PIN LOCK HOLE</b> Titanium	3.2 mm 38 g	<b>TAH T</b>
 <b>THREADED ADAPTER W/PIN LOCK HOLE</b> Stainless Steel	3.2 mm 66.4 g	<b>TAH S</b>
 <b>THREADED RECEIVER</b> Titanium	13.2 mm 36.8 g	<b>TR T</b>
 <b>THREADED RECEIVER</b> Stainless Steel	13.2 mm 64.5 g	<b>TR S</b>
 <b>3-PRONG ANCHOR, ROTATABLE</b> Stainless Steel	9.9 mm 122 g	<b>3PA R</b>
 <b>4-PRONG ANCHOR, ROTATABLE</b> Stainless Steel	16.3 mm 156 g	<b>4PA R</b>




COMPONENT	BUILD HT WGHT	PART NUMBER
 <b>4-HOLE PLATE, SQUARE, THREADED</b> Aluminum	6.4 mm 36 g	<b>4HLP ST A</b>
 <b>4-HOLE PLATE, ROUND, THREADED, TIE GROOVE</b> Aluminum	9.7 mm 59 g	<b>4HLP RTG A</b>
 <b>4-HOLE PLATE, ROUND, TIE GROOVE</b> Aluminum	9.7 mm 62 g	<b>4HLP RG A</b>
 <b>4-HOLE PLATE, ROUND, THREADED</b> Aluminum	6.4 mm 54 g	<b>4HLP RT A</b>
 <b>4-HOLE PLATE, SQUARE</b> Aluminum	6.4 mm 40 g	<b>4HLP A</b>
<b>WEIGHT LIMIT: 400 lbs (181 kg)</b>		
 <b>4-HOLE PLATE, ROUND, HEAVY DUTY</b> Titanium	6.4 mm 36 g	<b>4HLP RHD T</b>
 <b>4-HOLE PLATE, ROUND, HEAVY DUTY</b> Stainless Steel	6.4 mm 68 g	<b>4HLP RHD S</b>
<b>WEIGHT LIMIT: 165 lbs (75 kg)</b>		
 <b>4-HOLE PLATE, SQUARE, PEDIATRIC</b> Aluminum	6.4 mm 31 g	<b>4HLP AP</b>





# LAMINATION

*kits*

**MATERIAL**  
MULTI

**WEIGHT LIMIT**  
300 lbs (136 kg)

KIT	BUILD HT WGHT	PART NUMBER
 <b>3 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE</b> Titanium	13.1 mm  159.5 g	<b>3PA AHR T</b>
 <b>3 PRONG ANCHOR W/ RECEIVER, ROTATABLE</b> Titanium	23.1 mm  158.4 g	<b>3PA RR T</b>
 <b>4 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE</b> Titanium	19.4 mm  194 g	<b>4PA AHR T</b>
 <b>3 PRONG ANCHOR W/ RECEIVER, ROTATABLE</b> Titanium	29.5 mm  193 g	<b>4PA RR T</b>




KIT	BUILD HT WGHT	PART NUMBER
 <b>3 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE</b> Stainless Steel	13.1 mm  188 g	<b>3PA AHR S</b>
 <b>3 PRONG ANCHOR W/ RECEIVER, ROTATABLE</b> Stainless Steel	23.1 mm  186.1 g	<b>3PA RR S</b>
 <b>4 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE</b> Stainless Steel	19.4 mm  222 g	<b>4PA AHR S</b>
 <b>4 PRONG ANCHOR W/ RECEIVER, ROTATABLE</b> Stainless Steel	29.5 mm  220.5 g	<b>4PA RR S</b>




# LAMINATION

*additional parts*

**MATERIAL**  
MULTI

**WEIGHT LIMIT**  
MULTI

COMPONENT	WEIGHT	PART NUMBER
 <b>LAMINATION PLUG FOR ANCHOR</b> Nylon	24.6 g	<b>LP A</b>
 <b>LAMINATION CAP FOR 4-EAR SOCKET ADAPTER</b> Polyurethane	17.7 g	<b>LC</b>
 <b>LAMINATION PLUG FOR 3/4\"</b> 16 THREADED HOLE Stainless Steel	35 g	<b>LP TH</b>



COMPONENT	CONTENTS	PART NUMBER
 <b>LAMINATION PLATE FASTENER KIT</b> (4pc) M6 x 6 mm	Four fasteners (6 mm)	<b>LPFK 6</b>
 <b>LAMINATION PLATE FASTENER KIT, PEDIATRIC</b> (4pc) M5 x 6mm	Four fasteners (6 mm)	<b>LPFK 6P</b>
 <b>LAMINATION PLATE FASTENER KIT</b> (4pc) M6 x 10 mm	Four fasteners (10 mm)	<b>LPFK 10</b>

# FLEXION PLATES

**MATERIAL**  
MULTI





**WEIGHT LIMIT**  
300 lbs (136 kg)

KIT	BUILD HT WGHT	PART NUMBER
 <b>4-HOLE THREADED, 10°, W/ ADAPTER W/ PIN- LOCK HOLE, ROTATABLE</b> Titanium & Aluminum	12.9 mm   168 g	<b>4HTP 10 AHR T</b>
 <b>4-HOLE THREADED, 10°, W/ RECEIVER, ROTATABLE</b> Titanium & Aluminum	22.9 mm   166.8 g	<b>4HTP 10 RR T</b>
 <b>4-HOLE THREADED, 10°</b> Aluminum	9.7 mm   130 g	<b>4HTP 10 A</b>

KIT	BUILD HT WGHT	PART NUMBER
 <b>4-HOLE THREADED, 10°, W/ ADAPTER W/ PIN- LOCK HOLE, ROTATABLE</b> Stainless Steel & Aluminum	12.9 mm   196.3 g	<b>4HTP 10 AHR S</b>
 <b>4-HOLE THREADED, 10°, W/ RECEIVER, ROTATABLE</b> Stainless Steel & Aluminum	22.9 mm   194.5 g	<b>4HTP 10 RR S</b>

## ADDITIONAL HARDWARE

### STANDARD

HARDWARE	CONTENTS	PART NUMBER
 <b>4-HOLE FASTENER KIT</b> (4pc) M6 x 14 mm, for Non-Rotatable Adapters	Four fasteners (14 mm)	<b>4HFK 14</b>
 <b>4-HOLE FASTENER KIT</b> (4pc) M6 x 16 mm, for Rotatable Adapters	Four fasteners (16 mm)	<b>4HFK 16</b>
 <b>SUPER SET SCREWS</b> (8pc) M8 x 14 mm	Eight Super Set Screws	<b>CP SSET</b>
 <b>TUBE CLAMP SCREW</b> M5 x 14 mm	One tube clamp screw	<b>CP TCS</b>

HARDWARE	CONTENTS	PART NUMBER
 <b>4-HOLE FASTENER KIT</b> (4pc) M6 x 22 mm	Four fasteners (22 mm)	<b>4HFK 22</b>
 <b>FLEXION PLATE CLAMP SCREW</b> M6 x 30 mm	One flexion plate clamp screw	<b>CP FPSC</b>

# COMPONENT ORDER FORM

## BILLING AND SHIPPING INFO

TO AVOID SHIPPING DELAYS, PLEASE COMPLETE EVERY SECTION

CUSTOMER #

P.O. #

DATE

### BILLING

FACILITY/ATTN

ADDRESS

CITY

STATE(PROV)/ZIP

PHONE

FAX

**CARRIER\*** ☐ UPS ☐ OTHER

\*ORDERS SHIP UPS GROUND ON THE FOLLOWING DAY UNLESS SPECIFIED

### SHIPPING (Leave blank if same as billing)

FACILITY/ATTN

ADDRESS

CITY

STATE(PROV)/ZIP

PHONE

FAX

**DATE REQUIRED**

TIME

☐ AM ☐ PM

PART NUMBER	QTY

PART NUMBER	QTY

## NOTES

000 FRM COF NA 170808

EMAIL COMPLETED ORDERS TO: [orders@college-park.com](mailto:orders@college-park.com)

### RETURN POLICY

Unopened components may be returned within 90 days of purchase. A Return Authorization is required for all returns. Please feel free to contact Customer Service with any questions or concerns. Details and warranty information is available at: [www.college-park.com/warranty](http://www.college-park.com/warranty)

[www.college-park.com](http://www.college-park.com)



**college park**

27955 College Park Dr, Warren, MI 48088, USA  
T 586.294.7950 | F 586.294.0067 | [college-park.com](http://college-park.com)



**college park**

800.728.7950 | 586.294.7950 | [www.college-park.com](http://www.college-park.com)