

Indications For Use

I. Metatarsus Adductus/Varus (MTA)

The device is used for the primary correction of MTA, in place of serial castings and corrective shoes, or for the maintenance of correction after castings. Either the MTA or the TEV-CRB series may be used, depending on the preference of the physician. The MTA series is available from birth to 8 months of age. The TEV-CRB series is available from birth to 4 years of age. While the lower component alone is effective in most cases, some physicians may prefer to add the upper component for better control of the leg (see [III](#)).

II. Congenital Clubfoot

The device is used for the maintenance of correction after castings or surgery for clubfoot, by controlling forefoot adductovarus, hindfoot varus, and ankle equinus. The TEV-CRB series should be used. The upper component is optional, depending on the severity of the residual ITT (see [III](#)).

III. Internal Tibial Torsion

The device is used for the primary treatment of persistent tibial torsion in children up to age 4, in place of the D-B bar. For ITT, both upper and lower components are used together.