


1323 Knee – Polycentric 5 Bar



The 1323, is designed to be a very stable 5 Bar Knee with a large stance flexion angle (up to 12°). Weighing in at less than 1KG and having extensive alignment adjusting abilities give the 1323knee a smooth and natural gait pattern.



Lotus Adapter



The 1323 knee joint is a polycentric knee with stance flexion control for maximum safety during stance and is best suited to be used for low to high K3 level. Practitioner adjustable from extreme stability to high stance flexion. The 1323 offers a pneumatic cylinder with independent flexion extension adjustment for swing phase. Can be used for either above knee or knee disarticulation.

The frame construction is made of superlight aluminum alloy and the 5-bar linkage is made of aircraft alloy. All axes have built-in bearing which allow for ultra smooth walking movement. The proximal attachment has A-P slide and Rotation adjustment. Optional adapters are pyramid adapter (for AK) and lotus adapter (for KD).

P/N	1323AP Pyramid Adapter	1323AL Lotus Adapter
Weight Limit	125KG / 275lbs	125KG / 275lbs
Build Height	177.5mm	183.8mm
Knee Center to Top	24mm	30.3mm
Knee Flexion	135 Degrees	135 Degrees
Product Weight	976g	998g
Material	Aluminum	Aluminum
Warranty	2 year	2 year


REPLACEMENT PARTS

P/N	Description
PA13	Pyramid Adapter for: 1321, 1322, 1323, 1324 and 1326
LA13	Lotus Adapter for: 1321, 1322, 1323, 1324 and 1326



PDAC approved HCPCS code(s)
• L5840, L5845, L5850

To download PDAC Letter please visit our website at www.stngco.com



Stance Flexion Control (SFC)

SFC applies to L1324, 1324 & 1323 knees



ST&G USA CORPORATION
2781 Saturn St., Suite K
Brea, CA 92821

TELEPHONE: (714) 524-0663

FAX: (714) 364-8113

E-MAIL: info@stngco.com

www.stngco.com