


Genu Dynastab®
PDAC APPROVED L1810

Recommended for: Functional treatment of mild sprains and the resumption of sports activities after a sprain (mild, moderate, or severe)

- Dual hinge design provides stability
- Semi-rigid strapping gives extra support and prevents migration
- Patented pull tabs for easy application and removal
- Universal left or right


	Size	1	2	3	4	5	6
 Knee CIRCUM		12 ½ – 13 ¾"	13 ¾ – 15"	15 – 16 ½"	16 ½ – 18"	18 – 19 ½"	19 ½ – 21"


Genu Ligaflex® ROM
PDAC APPROVED L1832/L1833

Recommended for: Mild to moderate ligament injuries and instabilities / early stage OA / meniscus injuries and post-operative meniscal repairs.

Features

- 16" length with anterior closure
- Medical grade compression
- Easy to use Velcro tabs
- Low profile tension buckles
- TM5+ Hinge that replicates the natural movement of the knee
- Patented integrated anti-slip silicon knitting
- Tool-less hinge with easy to apply flexion and extension stops to restrict range of motion
- Contoured aluminum uprights to improve fit and comfort
- Universal left or right

	Size	1	2	3	4	5	6
 Knee CIRCUM		12 ½ – 13 ¾"	13 ¾ – 15"	15 – 16 ½"	16 ½ – 18"	18 – 19 ½"	19 ½ – 21"



Thuasne knit technology
Comfort, compression and breathability



Genu Ligaflex®
PDAC APPROVED L1820

Recommended for: Resumption of non-contact sports and regular daily activities after mild or moderate knee injuries.

Features

- Anterior wrap closure model: ultra easy to fit with speed clips
- Medical grade compression
- Hyper-flexible knitted areas in the popliteal crease and proximal edge
- Semi-rounded straps for improved comfort
- Thuasne's patented TM5 Hinge technology
- Dual tab calf strap improves suspension
- Integrated anti-slip silicone in the thigh section
- Aluminum uprights provide lateral stability
- Double adjustable hook and loop tabs help optimize placement of uprights

Genu Ligaflex Anterior Closure

	Size	1	2	3	3+	4	4+	5	5+	6	6+
 CIRCUM 5" Above Knee		14 ¾ – 17 ¼"	16 – 18 ¼"	17 ¼ – 19 ¾"	19 – 21 ¼"	19 – 21 ¼"	20 ½ – 22 ¾"	20 ½ – 22 ¾"	22 – 24 ½"	22 – 24 ½"	23 ½ – 26"
 CIRCUM 5 ½" Below Knee		11 ¾ – 13 ½"	13 – 14 ½"	14 ½ – 16"	14 ½ – 16"	16 – 17 ½"	16 – 17 ½"	17 ½ – 19"	17 ½ – 19"	19 – 20 ¾"	19 – 20 ¾"





A



B



C




D

Patella Reliever (A)
PDAC APPROVED L1810/L1812

Recommended for: Patello-femoral pain syndrome / lateral patellar subluxation / lateral patellar dislocations / patello-femoral arthritis / lateral patellar compression syndrome


- Medical grade compression
- Crescent shaped patella insert with lateral tracking system
- Distal strap for improved suspension
- Patented integrated antislip silicone knitting in the upper thigh
- Available in left and right

	Size	1	2	3	4	5	6
 Knee CIRCUM		12 ¼ – 13 ½"	13 ¾ – 14 ½"	14 ½ – 15 ¾"	15 ¾ – 17"	17 – 18"	18 – 19 ¼"

Ligastrap Genu (C)
PDAC APPROVED A4466

Recommended for: Maintenance, compression and support of the knee after mild or moderate knee strains or sprains


- Medical grade compression
- Numbered figure eight strapping delivers extra stabilization
- Anatomically shaped elastic knit for optimal thigh, calf, and popliteal comfort
- Patented pull tabs for easy application and removal
- Universal left or right

	Size	1	2	3	4	5	6
 Knee CIRCUM		12 ½ – 13 ¾"	13 ¾ – 15"	15 – 16 ½"	16 ½ – 18"	18 – 19 ½"	19 ½ – 21"

Silistab Genu (B)
PDAC APPROVED A4467

Recommended for: Maintenance / compression and support of the knee after mild or moderate knee strains or sprains


- Medical grade compression
- Patented silicone patella insert provides kneecap support and protects the patella tendon
- Breathable cotton fabric integrated posteriorly
- Flexible lateral stays give additional reinforcement
- Anatomically shaped elastic knit for optimal thigh, calf, and popliteal comfort
- Lateral silicone anti-slip system for better hold
- Universal left or right

	Size	1	2	3	4	5	6
 Knee CIRCUM		12 ½ – 13 ¾"	13 ¾ – 15"	15 – 16 ½"	16 ½ – 18"	18 – 19 ½"	19 ½ – 21"

Genuextrem (D)
PDAC APPROVED A4467

Recommended for: Maintenance, compression and support of the knee after mild or moderate knee strains or sprains.

- Medical grade compression
- Flexible lateral stays for additional reinforcement
- Anatomically shaped elastic knit for optimal thigh, calf, and popliteal comfort
- Lateral silicone anti-slip system for better hold
- Patented pull tabs for easy application and removal
- Universal left or right

	Size	1	2	3	4	5	6
 Knee CIRCUM		12 ½ – 13 ¾"	13 ¾ – 15"	15 – 16 ½"	16 ½ – 18"	18 – 19 ½"	19 ½ – 21"