


## Knee Joints



Lotus Adapter



Lanyard for Manual Locking Knee



### L1321

#### POLYCENTRIC 4-BAR MANUAL LOCKING KNEE

- Manual locking polycentric knee
- Lanyard lock mechanism can be disabled at a later time to convert to a non-locking free swing knee
- Provides a high degree of safety
- Frame construction is made of superlight aluminum alloy
- 4-bar linkage is made of aircraft alloy (AA7075)
- Proximal attachment has 2mm of A/P, M/L slide and rotation adjustment
- Optional Adapters: Pyramid (AK) and Lotus (KD)

#### PDAC approved HCPCS code(s)

- **L5850 + L5611 + L5925**
- To download PDAC Letter please visit our website at [www.stngco.com](http://www.stngco.com)

BILLING LEVEL **K1** **K2** **K3** **K4**

P/N	L1321AP Pyramid Adapter	L1321AL Lotus Adapter
<b>Weight Limit</b>	125kg / 275lbs	125kg / 275lbs
<b>Build Height</b>	108mm	114,3mm
<b>Knee Center to Mid Dome</b>	23mm	29,3mm
<b>Product Weight</b>	651g	672g
<b>Material</b>	Aluminum	Aluminum
<b>Warranty</b>	2 years	2 years

#### REPLACEMENT PARTS

P/N	Description
PA13	Pyramid Adapter for : L1321, 1321,1322,1323,1324, L1324 and 1326
LA13	Lotus Adapter for : L1321, 1321,1322,1323,1324, L1324 and 1326
LY13	Lanyard for Manual Locking Knee

Torque Value	Torque (Nm)	Hex (mm)
M6 Cap Screw (Tube Clamp Torque Setting)	12 Nm	5mm wrench
M10 Bolt (Pyramid Center Bolt)	18 Nm	8mm wrench