

**Newport® jr Pediatric Hip Orthosis** The Newport jr Hip System with the new pj—Virtual Pediatric Joint enhances hip stability by supporting and positioning the hips in optimal alignment while allowing controlled range of motion. This pediatric hip orthosis can be used post-operatively to maintain the corrected surgical position, or can be provided to maximize function during activities of daily living. Improved alignment through the hip encourages appropriate acetabular and femoral modeling, and reduces the tendency for subluxation and dislocation.

#### Indications


- Post-operative stabilization of the pediatric hip
- Patients with subluxed hips who require better alignment to relieve pain
- Children needing improved sitting and standing balance

#### Features

- Universal pelvic module, bilateral thigh cuffs, comfortable lumbar pad
- Anatomical design captures and stabilizes the pelvis
- Washable Coolfoam™ liners
- Available in five sizes

#### Virtual® Pediatric Joint (pj)

- The smallest and lightest hip joint available in the industry
- Adjustable control of flexion, extension, abduction and adduction
- Curved proximal joint attachments provide clearance over the incision and trochanter
- Joints with straight bars available in short, long and unslotted configurations for sale as components



Size	Model #	Waist Circ.	Hip Circ.	Mid Thigh Circ.	Trochanter to Knee Center Length	pj Joint Distal Bar Size	pj Joint Model #
X-Small	<b>3919</b>	17 <sup>7</sup> / <sub>8</sub> –19"	18–20"	9 <sup>7</sup> / <sub>8</sub> –10 <sup>7</sup> / <sub>8</sub> "	7 <sup>1</sup> / <sub>4</sub> –8 <sup>1</sup> / <sub>2</sub> "	Ex-short	<b>4310</b>
Small	<b>3920</b>	18 <sup>3</sup> / <sub>4</sub> –20 <sup>1</sup> / <sub>4</sub> "	20–22"	10 <sup>5</sup> / <sub>8</sub> –12 <sup>3</sup> / <sub>16</sub> "	8 <sup>1</sup> / <sub>8</sub> –10 <sup>1</sup> / <sub>8</sub> "	Short	<b>4311</b>
Medium	<b>3921</b>	20–22"	22–25"	11 <sup>3</sup> / <sub>16</sub> –13 <sup>3</sup> / <sub>4</sub> "	9 <sup>1</sup> / <sub>2</sub> –11 <sup>7</sup> / <sub>8</sub> "	Short	<b>4311</b>
Large	<b>3922</b>	21 <sup>3</sup> / <sub>4</sub> –23 <sup>3</sup> / <sub>4</sub> "	25–27"	12 <sup>3</sup> / <sub>4</sub> –15 <sup>3</sup> / <sub>16</sub> "	11 <sup>1</sup> / <sub>2</sub> –14"	Long	<b>4312</b>
X-Large	<b>3923</b>	23 <sup>1</sup> / <sub>2</sub> –25 <sup>1</sup> / <sub>2</sub> "	27–30"	14–16 <sup>7</sup> / <sub>8</sub> "	13 <sup>7</sup> / <sub>8</sub> –16 <sup>3</sup> / <sub>8</sub> "	Long	<b>4312</b>

Add **RL** suffix to part numbers for replacement liner sets, (e.g., **3920RL** is a complete liner set for the small Newport jr and **3940RL** is a liner set for the small left Newport jr Thigh) Replacement liners available. Available unilateral too.

#### Custom Hinge Components (Straight Bars)

Size	Model #	Joint Center To Proximal Bar	Joint Center To End Of Distal Bar
Short	<b>4315</b>	3 <sup>13</sup> / <sub>16</sub> "	7 <sup>13</sup> / <sub>16</sub> "
Long	<b>4316</b>	3 <sup>13</sup> / <sub>16</sub> "	10 <sup>15</sup> / <sub>16</sub> "
Long—unslotted	<b>4317</b>	5 <sup>5</sup> / <sub>16</sub> "	10 <sup>15</sup> / <sub>16</sub> "

#### Pelvic and Thigh Components

Size	Pelvic Assembly Model #	Thigh Assembly Model #	
		Left	Right
X-Small	<b>3929</b>	<b>3939</b>	<b>3944</b>
Small	<b>3930</b>	<b>3940</b>	<b>3945</b>
Medium	<b>3931</b>	<b>3941</b>	<b>3946</b>
Large	<b>3932</b>	<b>3942</b>	<b>3947</b>
X-Large	<b>3933</b>	<b>3943</b>	<b>3948</b>

#### Optional Adjustable Abduction Bar

Size	Attached	Kit Only
Short (For XS and S Pelvic Sections—6.35" to 9.85" adjustment range)	<b>3950</b>	<b>3950.01</b>
Long (For M, L and XL Pelvic Sections—9.675" to 17" adjustment range)	<b>3951</b>	<b>3951.01</b>

