



**PACELINE®**  
ADVANCED MEDICAL SOLUTIONS

# Closed-End PVA Bags

A Superior Fabrication Tool



---

Closed-End PVA Bags  
help you create a perfect  
finish with ease, in less than  
half the time, even on the  
most difficult models.

Increase your productivity by dramatically  
reducing time spent preparing a model  
for lamination.

Closed-End PVA Bags are available through one of your preferred distributors.

# Closed-End PVA Bags

A Superior Fabrication Tool

The Closed-End PVA Bag is a superior fabrication tool that allows you to eliminate fabrication steps...

Saving You Time & Money!

- Single Step Cast Preparation
- Saves Time and Money
- Works on a Wide Variety of Shapes and Sizes

## SIZING AND ORDERING INFORMATION STOCK BAGS (5 PER PACKAGE)

CLOSED-END PVA SIZING CHART		
Closed-End PVA Size	Proximal Measurement	Distal Measurement
4 - 0	6" - 9"	2"
6 - 1	9" - 12"	3"
6 - 3		4" - 6"
6 - 5		7" - 10"
8 - 3	12" - 16"	4" - 6"
8 - 5		7" - 10"
10 - 3	16" - 20"	4" - 6"
10 - 5		7" - 10"
10 - 7		11" - 14"
12 - 3	20" - 24"	4" - 6"
12 - 5		7" - 10"
12 - 7		11" - 14"
14 - 3	24" - 28"	4" - 6"
14 - 5		7" - 10"
16 - 5	28" - 32"	7" - 10"
16 - 7		11" - 14"

EXAMPLE: If your proximal size is 15", and the distal size is 9", then the correct Closed-End PVA size would be an 8-5.



Scan to see it in use .

## Starter Kit



Contains 5 of the most commonly used sizes

