a curocton	HEKO FORM					SUBMIT FORMS Fax 866.868.1890 E-mail orders@surestep.net Customer Service 877.462.0711			PATTERNS	
<i>∴</i> surestep	Date / /								& STRAPS CLICK OR SCAN	
PATIENT	STY	YLE						NOTES		
Patient ID			Weight: Patients	un to 100 lbs						
Age (in years)		Custom	Weight. Futients	up to 100 ibs						
□ Male □ Female		Custom	□ Polypropylene	(white)	Pattern:			S		
Weight Height			□ LDPE		Color: Black	□ Purple		NOTES		
Dx		□ Prefab	Weight: Patients	ht: Patients up to 50 lbs						
□ Right □ Left □ Bilateral			Only available in	y available in LDPE Color: □ Black		T.D. ale				
COMPANY			•	LDFE	Color: Black	□ Purple				
Name	M	EASUREM	ENTS							
Contact Name		Circumfe	ronco					$ \setminus \setminus $	\	P
Phone Fax			rence						<u> </u>	l Õ
E-mail	_	Distance M/L			\			$\prec \sim \sim$		
BILLING / SHIPPING		A/P		<u> </u>	_, l	•			<u> </u>	PROXIMAL TO KNEE
OrtoPed ULC 373 McCaffrey • Saint-Laurent, QC H4T 1Z7	ALL MEASUREMENTS					(/ 6	To
									Ş	
	SHOULD BE TAKEN IN WEIGHT BEARING				⊣	,	(\	—— 4"	血血
Shipping address same as billing? ☐ Yes ☐ No		WEIGHT E	SEARING	KNFFC	ENTER TO	[
Address City			TROCHAN		\		<u></u>	2"	CENTER	
	DF:	SIRED FINI								
ST/Prov Zip		NGTH				$H \times H$				~
PO#	Inch	nes proximal	to knee center:	<u> </u>		サムル	*	*/	KNEI	=
Order confirmation: 🗆 Fax 🗆 E-mail		•						\rightarrow	CENTER***	
SHIPPING OPTIONS									,	
Need by Date	Inches distal to		knee center:			((0))		\wedge		
GroundFREE				<u> </u>	(2"	D	
24-Hour Rush Fabrication Additional Charges apply separately						1 11 11 1				
						\ \\\\ /				ISTAL
		* Measurements required for HEKO Prefab orders			<u> </u>	\ \\ 	(4"	. TC
CASTS					NTER TO	. \}}\\				<u>조</u>
Sending Casts: ☐ Yes ☐ No	** 8" only nece HEKO hinges *** If knee cen		essary for large	LATERAL	L MALLEOLUS					Ë
Return impressions pre-modifications \$20.00/pr.			tor			_ /v 4	(<u> </u>	<u>C</u>
· · · · ·	cir	cumference	is greater than	:]				Z
Send casts: 373 McCaffrey Saint-Laurent, QC H4T 1Z7	11" large HEKO hinges will be used (additional charges					/			\	TO KNEE CENTER
		apply).						/7	8" **	