

Pylons and Pylon Adapters

All pylons may be cut to length with a tube cutting guide or other equipment. Some pylons have a permanent adapter at one end while others come with 30 mm tube inserts. Tube inserts are included.

30 MM ALUMINUM TUBE WITH ALUMINUM ADAPTER



PRODUCT WEIGHT 298 g

PATIENT WEIGHT LIMIT* 250 lb (115 kg) for U.S. Level 1 & 2; 220 lb (100 kg) for U.S. Level 3

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
Aluminum Fixed Adapter with 30 mm Aluminum Tube, 400 mm	FND-131330

30 MM ALUMINUM TUBE WITH TITANIUM ADAPTER



PRODUCT WEIGHT 313 g

PATIENT WEIGHT LIMIT* 250 lb (115 kg) for U.S. Level 1 & 2; 220 lb (100 kg) for U.S. Level 3

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
Titanium Fixed Adapter with 30 mm Aluminum Tube, 400 mm	FND-221330

30 MM ALUMINUM TUBE



For use with transfemoral sockets.

PRODUCT WEIGHT 193 g

PATIENT WEIGHT LIMIT* 250 lb (115 kg) for U.S. Level 1 & 2; 220 lb (100 kg) for U.S. Level 3

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
30 mm Aluminum Tube, 400 mm (two inserts included)	FND-011330
30 mm Tube Insert	FND-000030

MAGNUM 30 MM COMPOSITE TUBE



PRODUCT WEIGHT 251 g

PATIENT WEIGHT LIMIT* 350 lb (160 kg) for U.S. Levels 1-3; 300 lb (136 kg) for U.S. Level 4

WARRANTY 12 months from date of patient fitting
24 months from date of patient fitting when purchased as part of a Magnum Heavy Duty Kit

ITEM NUMBERS

Product	Item #
Magnum 30 mm Composite Tube, 400 mm (two inserts included)	MHD-30CT
Magnum 30 mm Composite Tube Insert	MHD-30CTI
Magnum Tube Cutting Guide	700-TL003

*Body weight plus any loads normally or routinely carried cannot exceed this weight limit.

Pylon Adapters

Connecting a pylon to a knee or a receiver is the job of endoskeletal adapters. These adapters can be especially useful in hard-to-fit cases or when there is little clearance. The Magnum 30 mm adapter is actually taller to provide more clamping force for heavy duty applications. Set Screws are included with all pylon adapters.

ALUMINUM 30 MM PYLON ADAPTER



PRODUCT WEIGHT 86 g

PATIENT WEIGHT LIMIT* 250 lbs (115 kg)

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
Aluminum 30 mm Pylon Adapter	FND-130030

TITANIUM 30 MM PYLON ADAPTER



PRODUCT WEIGHT 100 g

PATIENT WEIGHT LIMIT* 250 lb (115 kg)

WARRANTY 12 months from date of patient fitting

ITEM NUMBERS

Product	Item #
Titanium 30 mm Pylon Adapter	FND-220030

MAGNUM TITANIUM 30 MM ADAPTER



PRODUCT WEIGHT 110 g

PATIENT WEIGHT LIMIT* 350 lb (160 kg) for U.S. Levels 1-3; 300 lb (136 kg) for U.S. Level 4

WARRANTY 12 months from date of patient fitting
24 months from date of patient fitting when purchased as part of a Magnum Heavy Duty Kit

ITEM NUMBERS

Product	Item #
Magnum Titanium 30 mm Adapter	MHD-002

*Body weight plus any loads normally or routinely carried cannot exceed this weight limit.