



LaunchPad [launchpad-op.com](http://launchpad-op.com)

# X-tension

**X-tension** is first to market with the development of post fitting/fabrication dorsiflexion limiter for custom AFOs. It is designed to prevent drop off at terminal stance by resisting or stopping tibial progression over the foot.

X-tension bands generate plantar flexion assistance resulting in; increased third rocker power, contralateral step length, and single limb stance stability.

X-tension bands come in dynamic (flexible) and static (rigidly locks the AFO) options in both pediatric and adult sizes.


As part of LaunchPads TRIO innovations, X-tension works in combinations with SNAPstop and PIVOT ankle joints to achieve desired shank-to-floor alignments and tuning optimization.

Consider treating those sagittal plane compensations with knee extension, and flexion deficits, by increasing single limb stance stability and overall ambulatory confidence in your patient.

**LCODE recommendation:** L2220 x 2


Dorsiflexion and Plantar flexion assist/resist

## SIZES (in centimeters)



## DYNAMIC SIZES (in inches)

1/8", 5/32", 3/16", 1/4"



STATIC\*

DYNAMIC



\*—Static is available in one thickness



## HOW TO ORDER

**Cascade**  
 The Largest Independent  
 O&P Distributor

**Phone:** 800-888-0865

**Website:** [www.cascade-usa.com](http://www.cascade-usa.com)

**Email:** [orders@cascade-usa.com](mailto:orders@cascade-usa.com)