

# ALPS™ Prosthetic Ointment

P0840

## Prosthetic Ointment Soothes and Protects



ALPS Prosthetic Ointment is a vitamin A and D formulation that soothes and protects the skin. It may be applied nightly after the prosthesis has been removed, or during the day to the residual limb under a silicone liner. In order to reduce shear forces ointment can be applied at the proximal edge of the liner, fibular heads, patella, or condyles if necessary.

ALPS Prosthetic Ointment is available in a 4 oz. tube.

Recommended for use with:

- ✗ Gel (thermoplastic) - All ALPS Gel Liners\*
- ✓ Silicone - CLEARPRO
- ✓ Urethane

\*May deteriorate gel over time

# ALPS™ Prosthetic Anti-Perspirant ATP008

## Prosthetic Anti-Perspirant



The first anti-perspirant truly designed for prosthetics. Created to help manage the continuing problem of perspiration, it has been formulated with twice the Aluminum Chlorohydrate found in over the counter anti-perspirants.

ALPS Prosthetic Anti-Perspirant is available in a 4 oz. bottle.

Recommended for use with:

- ✓ Gel (thermoplastic) - All ALPS Gel Liners
- ✓ Silicone - CLEARPRO
- ✓ Urethane

SKIN CARE

