

# Stack Finger Splints

## *Flesh Style*



\* Stack Finger Splints are ideally suited for the repair of ruptured tendons. They are designed to support the distal joint of the finger in extension while permitting unrestricted movement of the proximal interphalangeal joint. The Splint is normally used without padding and is fixed to the middle phalanx with a small strip of adhesive tape, thus automatically elevating the distal joint to its proper extended position.

Occasionally, the dorsum of the finger may be too tender to tolerate a normal size Stack Splint. A much larger size, padded on the inside, should be used until the swelling has subsided enough to replace it with a Splint of the proper size. Because they are made of lightweight, flesh-colored plastic, patients will find Stack Finger Splints comfortable to wear for long periods of time. \*H. Graham Stack, F.R.C.S.

## ***BIRD & CRONIN INC.***

*Designers and Manufacturers of Quality Medical Products*

1200 Trapp Road • Eagan, Minnesota USA 55121

651/683-1111 • Fax 651/688-9855 • 800/328-1095 • Fax 800/279-7934

e-mail [sales@birdcronin.com](mailto:sales@birdcronin.com) • [www.birdcronin.com](http://www.birdcronin.com)